


Thanks to navigation, X-Guide's **Navigated FastMap** export contains Scan Bodies in the same coordinate system as the Intraoral scan (IOS) used to digitally design the Lab's wax-up in the treatment plan.




IOS soft tissue scans



- This means **there is NO need to take Pre- and Post-surgical soft tissue IOS scans with fixation fiducial screws to align the photogrammetry capture to the Lab Digital Wax-up/Treatment Plan IOS, which is required for traditional (basic) photogrammetry units.**

- FastMap can capture from All on 3 to All on 8 MUA positions.




1

FastMap Navigated Photogrammetry: LAB Workflow At a Glance


Digital Wax-up is transferred throughout the workflow to maintain the same coordinate system.

- A.** Lab creates Digital Wax-up on Treatment Plan IOS.



This information is integrated into the X-Guide Plan.

- B.** The clinician performs the FastMap Navigated Photogrammetry Scan.




Immediately after Navigated Implant surgery. The Patient Tracker Pattern serves as the fixation fiducial screws, and locates the patient anatomy relative to the Treatment Plan and relative to the final implant placement.

- C.** The clinician Exports FastMap Open STL: Scan Bodies & Treatment Plan IOS in same coordinate system.



Export converts FastMap Scan Body to Elos IO 2C-A AND to Nobel White Caps.

- D.** Merge in Lab Design Software. No need for fiducial IOS scans.



Keep in mind, there is no need for IO scans with fixation fiducial screws as the X-Guide Tracker & Registration replaces it.

- E.** The Lab Finalizes Provisional Design for Fabrication.



2

FastMap Output data will come in a ZIP single file



[DOWNLOAD SAMPLE Navigated FASTMAP DATASET](#)

Uncompressed File Structure – Navigated FastMap:

Name	Type
Advanced	File folder
SupplementaryData	File folder
Export_Summary	Text Document
MAXILLA_Anatomy_IOS	STL File
MERGED_MAXILLA_Elos-IO2C-A-NobelBiocare-MultiUnitNPRP	STL File

Import these 2 files into design software. We auto convert the FastMap Scan Bodies into the Elos IO 2C-A geometry.

NOTE: Nobel White Cap & more Elos IO data is in the **ADVANCED** folder.

ADVANCED Folder:

Name	Type
_Elos-IO2C-A-NobelBiocare-MultiUnitNPRP	File folder
_Elos-IO2C-B-NobelBiocare-MultiUnitWP	File folder
_white-cap-NobelBiocare-MultiUnitNPRP	File folder
_white-cap-NobelBiocare-MultiUnitWP	File folder
_XNav-IO2C-A-NobelBiocare-MultiUnitNPRP	File folder
Diagnostics	File folder
MERGED_MAXILLA_Elos-IO2C-A-NobelBiocare-MultiUnitNPRP	STL File
MERGED_MAXILLA_Elos-IO2C-B-NobelBiocare-MultiUnitWP	STL File
MERGED_MAXILLA_white-cap-NobelBiocare-MultiUnitNPRP	STL File
MERGED_MAXILLA_white-cap-NobelBiocare-MultiUnitWP	STL File
MERGED_MAXILLA_XNav-IO2C-A-NobelBiocare-MultiUnitNPRP	STL File

FAQ's:

Do we need a post-surgical IOS soft tissue scan?

- It is not required for aligning the meshes in the design software, but if you would like a post-surgical soft tissue profile, the clinician can capture that IOS. The fixation fiducial screws are not required. If taking this scan, the clinician should use the Nobel Healing Cap Multi-unit (29064/5 pack or 31145/single pack) or Elos IO 2C-A during capture. (FastMap will export both geometries in the same coordinate system with the original pre-op treatment planning IOS used for digital wax-up.)



Nobel White Cap Healing Abutment
(For the Nobel Biocare workflow, please use their white cap healing abutment for the post surgical IOS soft tissue scan.)



or





Elos IO
2C-A

- Some labs choose to skip the post-surgical IOS scan (since it is not needed to merge datasets). Labs design the prosthesis assuming soft tissue is 2 or 3mm below MUA level. Lab and clinician must agree on design method.

FastMap® Navigated Photogrammetry for the Dental Lab

In Exocad, can I use DESS Screws with FastMap Data?

- Yes, this is how: Select the library parts from the Select Library drop down menu in the Wizard. Match to the scan abutment using the red marked spot.




FAQ's:

How do we access different geometry (shapes) libraries in Exocad Design Software?

Click [EXOCAD LINK](#) to download various packages for: **ELOS**, **Nobel Biocare**, **DESS**

ELOS

Download Packages: Elos IL IO exocad libraries (for non-US customers) and Elos IL FDA IO exocad libraries (for US customers). This package contains the Elos Accurate – Hybrid Base Non-Engaging library, which contains the Elos IO 2C-A STL geometry.




The screenshot shows the ELOS download interface. It features the ELOS logo on the left. On the right, there are three package cards, each with a title, update date, number of libraries included, and an 'Add to download list' button. The first card is for 'Elos Lab Analog exocad libraries' (30-08-2024, 4 libraries). The second card is for 'Elos IL IO exocad libraries' (06-03-2025, 102 libraries) and is highlighted with a red border. The third card is for 'Elos IL FDA IO exocad libraries' (12-02-2025, 41 libraries) and is also highlighted with a red border.

Package	Updated	Libraries Included
Package: Elos Lab Analog exocad libraries	30-08-2024	4
Package: Elos IL IO exocad libraries	06-03-2025	102
Package: Elos IL FDA IO exocad libraries	12-02-2025	41

Nobel Biocare

Download Packages: Nobel Biocare Non US IL ALL exocad libraries and Nobel Biocare IL FDA All exocad libraries. This package contains the Nobel Biocare – Temp Bridge library, which contains both the Elos IO 2C-A and Nobel whitecap (Healing Cap 31145) geometries with specific DESS screw sizes.




The screenshot shows the Nobel Biocare download interface. It features the Nobel Biocare logo on the left. On the right, there are five package cards, each with a title, update date, number of libraries included, and an 'Add to download list' button. The first card is for 'NobelBiocare Non US IL All exocad libraries' (12-12-2024, 11 libraries) and is highlighted with a red border. The second card is for 'NobelBiocare NobelPearl Lab Analog exocad libraries' (07-11-2024, 1 library). The third card is for 'NobelBiocare NobelPearl IL exocad libraries' (02-08-2024, 1 library). The fourth card is for 'NobelBiocare IL FDA All exocad libraries' (27-11-2024, 10 libraries) and is highlighted with a red border. The fifth card is for 'NobelBiocare FDA Lab Analog exocad libraries' (07-11-2024, 8 libraries).

Package	Updated	Libraries Included
Package: NobelBiocare Non US IL All exocad libraries	12-12-2024	11
Package: NobelBiocare NobelPearl Lab Analog exocad libraries	07-11-2024	1
Package: NobelBiocare NobelPearl IL exocad libraries	02-08-2024	1
Package: NobelBiocare IL FDA All exocad libraries	27-11-2024	10
Package: NobelBiocare FDA Lab Analog exocad libraries	07-11-2024	8

DESS

Download Packages: DESS FDA IO IL exocad libraries (US customers) and DESS Angled Base All IO IL exocad libraries and DESS BASE All IO IL exocad libraries (non-US customers).



The screenshot shows the DESS download interface. It features the DESS logo on the left. On the right, there are five package cards, each with a title, update date, number of libraries included, and an 'Add to download list' button. The first card is for 'DESS PIC All IO exocad libraries' (22-10-2024, 4 libraries). The second card is for 'DESS Lab Analog exocad libraries' (28-01-2025, 32 libraries). The third card is for 'DESS FDA IO IL exocad libraries' (28-02-2025, 126 libraries) and is highlighted with a red border. The fourth card is for 'DESS Angled Base All IO IL exocad libraries' (04-03-2025, 84 libraries). The fifth card is for 'DESS BASE All IO IL exocad libraries' (28-02-2025, 51 libraries) and is also highlighted with a red border.

Package	Updated	Libraries Included
Package: DESS PIC All IO exocad libraries	22-10-2024	4
Package: DESS Lab Analog exocad libraries	28-01-2025	32
Package: DESS FDA IO IL exocad libraries	28-02-2025	126
Package: DESS Angled Base All IO IL exocad libraries	04-03-2025	84
Package: DESS BASE All IO IL exocad libraries	28-02-2025	51